MATCHPOST

Fiber Posts for Dental Restoration
Drills for Dental Restoration

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
NON-STERILE – STERILIZE BEFORE FIRST USE

INSTRUCTIONS FOR USE



PRODUCT DESCRIPTION

MATCHPOST devices consist of **MATCHPOST** posts and **MATCHPOST** drills.

MATCHPOST posts are dental posts used for coronal-radicular reconstruction of devitalized teeth. **MATCHPOST** posts are radiopaque translucent fiber posts. **MATCHPOST** posts are smooth-surfaced posts, with a predominantly parallel cylindrical shape, and a taper on the bottom quarter of the post. The length of **MATCHPOST** posts varies with their diametrical size, between 15-21 mm. The radio-opacity level is around 200% Al equivalent (ISO 4049).

MATCHPOST drills (STARTER DRILL/FINISHING DRILL) are intended for mechanical shaping of the root canal. The shape is chosen according to the form of the corresponding post: **MATCHPOST** drills exist in a similar range of diameters as the associated posts, and they are identified by a color code.

Posts and Drills sequential

INDICATIONS FOR USE

MATCHPOST posts have been designed to support and secure the coronal-radicular restoration in case of insufficient residual tooth structure (<4 mm).

MATCHPOST drills are intended to prepare the tooth for insertion of a post that will support and secure the coronal restoration, where insufficient residual tooth structure (<4 mm) exists.

These products are to be used only in a clinical or hospital environment, following Good Dental Practice, by qualified dental professionals such as general practitioners as well as Endo specialists (Endodontists) and Dental Assistants.

CONTRAINDICATIONS

MATCHPOST posts and associated drills are contraindicated if coronal tooth structure is less than 1.5 mm.

DELIVERY FORMS

The devices are supplied in:

- blister(s) containing 5 posts each;
- blister(s) containing 1 drill each.

The blisters may be supplied in different commercial presentations; see external packaging for content detail.

COMPOSITION

MATCHPOST posts are fiber-reinforced composite material, made of radiopaque fiber (80%) and epoxy resin matrix (20%).

MATCHPOST drills are made of stainless steel.

SAFETY NOTES

WARNINGS

Strictly follow these Instructions for Use for **MATCHPOST** posts and associated **MATCHPOST** drills to minimize the following residual risks to the device, the patient and/or the user:

- Root fracture, requiring tooth extraction
- Breakage (fracture) of device or instrument.
- Debonding, requiring new restoration
- Dentin damage due to heat generation in case of insufficient irrigation during drilling,
- Infection due to contamination or cross-contamination
- Swallowing or inhalation of device (e.g. post)
- Temporary discomfort; longer procedure.

Serious health consequences for the patient (e.g. biological hazards such as cytotoxicity, genotoxicity or carcinogenicity) have never been reported and are very unlikely but cannot be totally excluded.

PRECAUTIONS

There is no data for the safety and effectiveness of use of **MATCHPOST** posts and associated drills for pregnant women, breastfeeding women and children under 18 years old (for posts only). They are therefore not recommended for such uses.

Inspect the product before use and do not use the device(s) if the product (i.e. device, primary and/or secondary packaging and labeling) is damaged.

Do not use for indications different from the IFU's recommendations.

MATCHPOST posts:

- The post must be immersed in isopropanol for 1 minute then dried gently with air syringe for 15 seconds at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.
- Shortening the post with a diamond disc should be done outside the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post. Do not apply excessive pressure when seating the post in the root canal.
- The use of radiographs is recommended to assess root canal morphology, determine the post and drill size needed, and/or verify post position in the root canal.
- The use of a rubber dam is recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.
- For profoundly flared or ovoid canals, the use of accessory post* is indicated for best results.
- A selected re-access kit* can be used if needed.

MATCHPOST drills:

- Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull.
- Clean and sterilize the instrument as instructed prior to first use and between patients, and do not exceed the indicated maximum number of allowed uses.
- Verify proper and secure connection with contra-angle head before use and select appropriate rotation speed, as recommended. Do not apply excessive pressure on the drill.
- The STARTER DRILL may be used only to remove root canal filling partially, and the FINISHING DRILL may be used only to prepare the root canal.
- The use of a rubber dam is recommended.

ADVERSE REACTIONS / UNDESIRABLE SIDE EFFECTS

No known undesirable side effect.

STORAGE CONDITIONS

Store at ambient temperature.

STEP-BY-STEP INSTRUCTIONS

1. Posts are to be sterilized prior to use. Drills are to be sterilized prior to first use and between uses and are to be replaced after 15 uses. (see Processing section)
2. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
3. Select the correct **MATCHPOST** post size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
4. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the **MATCHPOST STARTER DRILL** selected (rotation speed 800-1200 RPM, Revolutions Per Minute) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable. Irrigate as necessary to prevent rotation-induced heating.
5. Shape the canal with the **MATCHPOST FINISHING DRILL** corresponding to the selected **MATCHPOST** post size, (rotation speed 1000-2000 RPM, Revolutions Per Minute), irrigating as necessary to prevent heating.
6. Try in the post, to check that the post is properly seated in the canal.
7. Shorten the post to its final length with a diamond disc, outside the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

8. Immerse or soak the post in isopropanol for 1 minute and then gently air dry for 15 seconds at least.
9. Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove excess water with paper points but leave the surface moist. Do not put the etchant in contact with the gingiva.
10. Apply 2 coats* of primer with a thin brush in the post space. Remove excess primer using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light cure the primer* inside the post space for 10 to 20 seconds.
11. Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
12. a. If using a dual-cure resin cement for both cementation and core build-up such as **CORECEM** (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light cure for 40 to 60 seconds. Continue applying the core build-up using a plastic form matrix*, or directly on the tooth.
12. b. If using separate cement and core build-up resins: Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove any excess and air-dry with a jet of air. Light cure for 20 seconds (if needed).
12. c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
12. c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove any excess and air dry with a jet of air. Light cure for 20 seconds (if needed).
13. Directly model the build-up using a composite core build-up material*. RTD recommends that the coronal end of the post be covered by core composite. In cases where a prosthetic crown is not indicated, a minimum of 1 mm of composite restorative material should cover the post, and that should be an appropriately shaded, polishable, aesthetic composite.

PROCESSING: HYGIENE, STERILIZATION

MATCHPOST posts:

- **MATCHPOST** post is a single use device. It must not be re-used to avoid the risk of contamination.
- Before insertion into the canal, the post must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g. EN ISO 11607-1),
 - ✓ Autoclave: type B complying with regulation standard in force (e.g. EN 13060),
 - ✓ Sterilization temperature: 134°C – sterilization time: 18 minutes, (e.g., “Prions” program)
 - ✓ 1 cycle only.
- **MATCHPOST** posts cannot be sterilized or disinfected in a thermoisinfectant.

MATCHPOST drills:

- **MATCHPOST** drills are supplied non-sterile. Disinfect and sterilize **MATCHPOST** drills before each use. Disinfect instruments with disinfecting or cleaning agents adapted to rotary instruments*. Do not disinfect instruments in a thermoisinfectant.
- The lifetime of use for the **MATCHPOST** drills is defined as 15 uses, following which, the practitioner or dental surgeon is informed that the drill should be replaced.
- **MATCHPOST** drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g. EN ISO 11607-1),
 - ✓ Autoclave: type B complying with regulation standard in force (e.g. EN 13060),
 - ✓ Sterilization temperature: 134°C – Sterilization: 18 minutes, (e.g., “Prions” program)
 - ✓ Store the sterilized components in a dry, dust-free place,
 - ✓ If the packaging integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

DISPOSAL

Dispose in a safe manner in accordance with local/national regulations for disposal of sharp and contaminated devices.

ADDITIONAL INFORMATION

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

WARRANTY

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

LIMITATION OF LIABILITY

Except where prohibited by law, RTD will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For Dental Use Only. For Professional Use Only.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Information available on Dentsply Sirona website: [Dentsply Sirona | USA \(https://www.dentsplysirona.com/en-us\)](https://www.dentsplysirona.com/en-us) and/or in the European database on medical devices (EUDAMED): <https://ec.europa.eu/tools/eudamed>, where it shall be linked to the Basic UDI-DI of **MATCHPOST posts**.

Basic UDI-DI numbers:


MATCHPOST posts: ++ERTD4800EUR2

MATCHPOST drills: ++ERTD4820EURG

*** Refer to the manufacturer's instructions for use.**

**** Please refer to your autoclave manufacturer for specific sterilization instructions if needed.**

TERMS AND SYMBOLS IDENTIFICATION (e.g. ACCORDING TO ISO 15223-1 AND ISO 21531)

	Medical Device	Posts and Drills
	Manufacturer	Posts and Drills
	Date of manufacture	Posts and Drills
	Do not re-use	Posts only
	Number of uses less one	Drills only
	Batch code	Posts and Drills

	Catalogue number	Posts and Drills
	Number of unit(s) / Content	Posts and Drills
	Do not use if package is damaged	Posts and Drills
	Consult Instructions for use	Posts and Drills
	Caution	Posts and Drills
	Professional Use	Posts and Drills
	Non-Sterile	Posts and Drills
	Sterilizable in a steam sterilizer (autoclave) at temperature specified (134°C)	Posts and Drills

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INSTRUCTIONS FOR USE


DESCRIPTION

MACRO-LOCK devices consist of **MACRO-LOCK** posts and **MACRO-LOCK** drills. **MACRO-LOCK** posts are radiopaque, translucent fiber post. **MACRO-LOCK** posts present a serrated surfaced post, apically tapered, with a length of 17.5 mm. The radio-opacity level is around 200% Al equivalent (ISO 4049).

MACRO-LOCK drills (STARTER/FINISHING) are intended for mechanical shaping and cleaning of the root canal. They exist in a similar range of diameters, as the associated post, and they are identified by a color code.

MACRO-LOCK drills can be associated with **MACRO-LOCK**, **MACRO-LOCK POST ILLUSION XRO** and **MACRO-LOCK OVAL** post devices.

	1	2	3	4
Ø APICAL TIP	0.80	0.80	1.00	1.00
Ø POST HEAD	1.35	1.47	1.67	1.83
STARTER DRILL	#1	#1	#2	#2
FINISHING DRILL	#1	#2	#3	#4

COMPOSITION

MACRO-LOCK post is made of 80% radiopaque quartz fibers and 20% epoxy resin matrix.

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INTENDED USE

MACRO-LOCK post has been designed to support and secure the coronal-radicular restoration of devitalized teeth.

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INDICATIONS FOR USE

MACRO-LOCK devices are indicated in case of insufficient residual tooth structure (<4 mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

MACRO-LOCK devices are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of **MACRO-LOCK** devices for pregnant women, breastfeeding women and children under 18 years old (posts only). It is not recommended for such use.

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MACRO-LOCK post:

- **MACRO-LOCK** post is single use device. It must not be re-used to avoid the risk of contamination.
- For profoundly flared or ovoid canals, the use of accessory post is indicated for best results.
- The post must be sterilized individually before insertion into the canal the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g EN ISO 11607-1),
 - ✓ Autoclave: type B complying with regulation standard in force (e.g EN 13060),
 - ✓ Sterilization temperature: 134°C- sterilization time: 18 minutes,
 - ✓ 1 cycle only.
- The post must be immersed in isopropanol for 1 min then dried gently with air syringe for 15s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- The use of rubber dam is recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.

MACRO-LOCK drills:

- **MACRO-LOCK** drills are supplied non-sterile. Disinfect and sterilize **MACRO-LOCK** drills before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfectant.
- Inspect **MACRO-LOCK** drills visually before each use: change to new drills when they are damaged, corroded or dull.
- The STARTER/FINISHING drill may be used only to remove root canal filling partially, and the FINISHING drills may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the **MACRO-LOCK** drills stated is 15 uses.
- **MACRO-LOCK** drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,
 - ✓ Autoclave: type B complying with regulation standard in force,
 - ✓ Sterilization temperature: 134°C – Sterilization: 18 minutes,
 - ✓ Store the sterilized components in a dry, dust-free place,
 - ✓ If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

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STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo or reamer drill torpan #0.70. Radiographic verification is recommended.
2. Select the correct **MACRO-LOCK** post size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for $\frac{1}{3}$ of the final post length. Remove the root filling with the **MACRO-LOCK** STARTER DRILL selected (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted to $\frac{2}{3}$ of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the **MACRO-LOCK** FINISHING DRILL corresponding to the selected **MACRO-LOCK** post size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal.
6. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
7. Immersed or soaked the post in isopropanol for 1 minute and then gently air dry for 15 s at least.
8. Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.
9. Apply 2 coats* of the primer with a thin brush in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth*.
- 11a. If using a dual cure resin cement for both cementation and core build up such as CORECEM (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up or directly on the tooth.
- 11b. If using separate cement and a core build up resins: Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
- 11c. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
- 12 Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the post be covered by core composite.

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SAFETY AND LIABILITY

WARRANTY

R.T.D. warrants this product will be free from defects in material and manufacture. R.T.D. makes no other warranties

including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and R.T.D.'s sole obligation shall be repair or replacement of the R.T.D. product.

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TERMS AND SYMBOLS IDENTIFICATION

	Manufacturer
	Do not re-use
	Batch code
	Catalogue number
	Consult instructions for use
	Do not use if package is damaged
	Caution
	Number of unit / Content

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INSTRUCTIONS FOR USE

DESCRIPTION

MACRO-LOCK POST ILLUSION X-RO are dental posts used for reconstruction of devitalized teeth. They are made of X-RO yellow translucent fibers bound in a colored resin Epoxy-vinyl ester matrix, which becomes translucent in the mouth. The post color will re-appear on command under cold air, water spray if re-treatment is required.

The MACRO-LOCK POST ILLUSION X-RO is a serrated-surfaced post, apically tapered with a length of 17.5 mm. The post has fluted slots around its axis for additional macro-retention to the adhesive used, by increasing the surface area available on the post for cement contact.

The radiopacity level of the MACRO-LOCK POST ILLUSION X-RO is 340 % Al equ. (ISO 4049)

The devices are available in 6 sizes and are used in association with a STARTER DRILL to perform a first hole and FINISHING DRILL is associated to each post size according to a color coding (see Figure 1). The drills associated to MACRO-LOCK POST ILLUSION X-RO are MACRO-LOCK drills.

	● 1	● 2	● 3	● 4	● 5	○ 6
Ø apical tip	0.80	0.80	1.00	1.00	1.30	1.30
Ø Post head	1.35	1.47	1.67	1.83	2.05	2.20
Starter drill	#1	#1	#2	#2	#2	#2
Finishing Drill	#1	#2	#3	#4	#5	#6


Figure 1 : MACRO-LOCK POST ILLUSION X-RO.

COMPOSITION

MACRO-LOCK POST ILLUSION X-RO posts are fiber-reinforced polymer, made of XRO fiber (80%) and epoxy resin matrix (20%) which contains colored pigment.

MACRO-LOCK drills are made of stainless steel.

INTENDED USE

MACRO-LOCK POST ILLUSION X-RO posts has been designed to support and secure the coronal restoration in case of insufficient residual tooth structure (<4 mm).

MACRO-LOCK drills are intended to prepare the tooth for insertion of a post that will support and secure the coronal restoration, where insufficient residual tooth structure (<4 mm) exists.

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INDICATIONS FOR USE

MACRO-LOCK POST ILLUSION X-RO posts and associated drills are indicated in case of insufficient residual tooth structure (<4mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

MACRO-LOCK POST ILLUSION X-RO posts and associated drills are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of **MACRO-LOCK POST ILLUSION X-RO** posts and associated drills for pregnant women, breastfeeding women and children under 18 years old (for posts only). It is not recommended for such use.

If the device is used in a way different than the one recommended in the IFU, loosening, breakage of the post or extraction of the tooth might occur.

MACRO-LOCK POST ILLUSION X-RO posts:

- The post must be immersed or soaked in isopropanol for 1 minute and dried gently with air syringe for 15 s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that. **The post can neither be sterilized nor disinfected in thermoisinfector and the use of autoclave is not recommended.**
- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- The use of rubber dam is recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.

MACRO-LOCK drills:

- **MACRO-LOCK POST ILLUSION X-RO** drills are supplied non-sterile. Disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfector.
- Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull.
- The STARTER DRILL may be used only to remove root canal filling partially, and the FINISHING DRILLS may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the **MACRO-LOCK** drills stated is 15 uses, following this the dental surgeon is informed that the drill should be replaced.
- Drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,

- ✓ Autoclave: type B complying with regulation standard in force,
- ✓ Sterilization temperature: 134°C - sterilization time: 18 minutes,
- ✓ Store the sterilized components in a dry, dust-free place,
- ✓ If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

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STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
2. Select the correct **MACRO-LOCK POST ILLUSION X-RO** size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the **MACRO-LOCK STARTER DRILL** selected (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the **FINISHING DRILL** corresponding to the selected **MACRO-LOCK POST ILLUSION X-RO** size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal.
6. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
7. Clean the post with alcohol.
8. **Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.**
9. Apply 2 coats* of the primer with a thin brush in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
11. a. If using a dual cure resin cement for both cementation and core build up such as **CORECEM** (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic form matrix, or directly on the tooth.
11. b. If using separate cement and a core build up resins : Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
11. c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
11. c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the post be covered by core composite.

STORAGE

Do not store the post in direct sun light or artificial light.

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SAFETY AND LIABILITY

WARRANTY

R.T.D. warrants this product will be free from defects in material and manufacture. R.T.D. makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R.T.D.'s sole obligation shall be repair or replacement of the R.T.D. product.

LIMITATION OF LIABILITY

Except where prohibited by law, R.T.D. will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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TERMS AND SYMBOLS IDENTIFICATION

	Manufacturer
	Batch code
	Catalogue number
	Number of unit / Content
	Keep away from sunlight

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CNPJ: 12.483.930/0001-22

www.dentsplysirona.com/pt-br - SAC: 0800 771 2226 / (11) 3046-2222

RTD fiber posts ILLUSION X-RO product family ANVISA: 80745409032

Universal/Starter Drills ANVISA: 80745400054

Finishing Drills: ANVISA: 80745400055

CAUTION: Make sure that this instruction for use corresponds to the product that you have purchased, as indicated on the packaging.

To obtain this instruction for use in a printed format free of charge, please call 0800 771 2226 (Brazil only) or send an email to atendimento@dentsplysirona.com

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INSTRUCTIONS FOR USE

DESCRIPTION

MACRO-LOCK OVAL are dental posts used for reconstruction of devitalized teeth. They are made of X-RO translucent fibers bound in a colored resin Epoxy-vinyl ester matrix, which becomes translucent in the mouth. The post color will reappear on command under cold air, water spray if re-treatment is required.

The **MACRO-LOCK OVAL** is a serrated surface post, apically tapered with a length of 18 mm. The post is round at the apical end and becomes oval on the coronal section.

The radiopacity level of the **MACRO-LOCK OVAL** is 340 % Al equivalent (ISO 4049).

The devices are available in 4 sizes and are used in association with **MACRO-LOCK STARTER/FINISHING** drills, as indicated on the picture here below.

MACRO-LOCK drills are associated with:

- **MACRO-LOCK** posts
- **MACRO-LOCK OVAL**
- **MACRO-LOCK POST ILLUSION X-RO**


Figure 1 : **MACRO-LOCK OVAL** posts and **MACRO-LOCK STARTER/FINISHING** drills.

COMPOSITION

MACRO-LOCK OVAL posts are fiber-reinforced polymer, made of XRO fiber (80%) and epoxy resin matrix (20%) which contains colored pigment.

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INTENDED USE

MACRO-LOCK OVAL posts have been designed to support and secure the coronal restoration in case of insufficient residual tooth structure (<4 mm).

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These products should only be used in hospitals, clinics or dental office with qualified dental personnel (dentist).

INDICATIONS FOR USE

MACRO-LOCK OVAL posts and associated drills are indicated in case of insufficient residual tooth structure (<4mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

MACRO-LOCK OVAL posts and associated drills are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of **MACRO-LOCK OVAL** and associated drills on pregnant women, breastfeeding women and children under 18 years old (for posts only). It is not recommended for such use.

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MACRO-LOCK OVAL posts:

- The post must be immersed or soaked in isopropanol for 1 minute and dried gently with air syringe for 15 s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that. **The post can neither be sterilized nor disinfected in thermoisinfector and the use of autoclave is not recommended.**
- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- The use of rubber dam is recommended and can momentarily prevent the post to change color.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.

MACRO-LOCK drills:

- **MACRO-LOCK** drills are supplied non-sterile. Disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfector.
- Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull.
- The **STARTER DRILL** may be used only to remove root canal filling partially, and the **FINISHING DRILLS** may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the **MACRO-LOCK** drills stated is **15 uses**, following this the dental surgeon is informed that the drill should be replaced.
- Drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,
 - ✓ Autoclave: type B complying with regulation standard in force,
 - ✓ Sterilization temperature: 134°C - sterilization time: 18 minutes,
 - ✓ Store the sterilized components in a dry, dust-free place,
 - ✓ If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

ADVERSE REACTIONS

- Debonding of the post

- Root fracture
- Fracture of the post

STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
2. Select the correct **MACRO-LOCK OVAL** size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the **MACRO-LOCK STARTER DRILL** selected (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the **MACRO-LOCK FINISHING DRILL** corresponding to the selected **MACRO-LOCK OVAL** size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal.
6. The contour of the post can be reduced slightly, if necessary, with a diamond bur. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
7. Clean the post with alcohol.
8. Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.
9. Apply 2 coats* of the primer with a thin brush in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
11. a. If using a dual cure resin cement for both cementation and core build up such as CORECEM (recommended technique): Mix the cement and apply it in the post space from bottom to the top using narrowest syringe tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic form matrix, or directly on the tooth.
11. b. If using separate cement and a core build up resins : Mix the cement* and apply it to the post and then in the post space using a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
11. c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
11. c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the post be covered by core composite.

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SAFETY AND LIABILITY

WARRANTY

R.T.D. warrants this product will be free from defects in material and manufacture. R.T.D. makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your

exclusively remedy and R.T.D.'s sole obligation shall be repair or replacement of the R.T.D. product.

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***Refer to the manufacturer's instructions for use.**

****Please refer to your autoclave manufacturer for specific sterilization instructions if needed.**

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38120 St Egrève - França

RTD-info@dentsplysirona.com; www.rtdental.com

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D.T. LIGHT-POST®

FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE

DESCRIPTION

D.T. LIGHT-POST® devices consist of D.T. LIGHT-POST® posts and D.T. LIGHT-POST® drills.



D.T. LIGHT-POST® posts are radiopaque, translucent fiber post. D.T. LIGHT-POST® posts present a smooth surfaced post, apically tapered with a double taper, and a length of 20mm. The radio-opacity level is around 200% Al equivalent (ISO 4049).

D.T. LIGHT-POST® drills (UNIVERSAL/FINISHER) are intended for mechanical shaping and cleaning of the root canal. They are made of stainless steel and the shape is chosen according to the form of the corresponding post. D.T. LIGHT-POST® drills exist in a similar range of diameters as the associated post, and they are identified by a color code.

D.T. LIGHT-POST drills are associated with:

- D.T. LIGHT-POST posts
- D.T. LIGHT-POST ILLUSION X-RO posts

	■ 0,5	■ 1	■ 2	■ 3
Ø apical tip	0.80	0.90	1.00	1.20
ØPost head	1.25	1.52	1.80	2.18
Double taper	.02/.04	.02/.06	.02/.08	.02/.10
Color coding	■ Black	■ Red	■ Yellow	■ Blue

COMPOSITION

D.T. LIGHT-POST® posts are fiber-reinforced polymer made of radiopaque quartz fiber (80%) and epoxy resin matrix (20%).

D.T. LIGHT-POST® drills are made of stainless steel.

INTENDED USE

D.T. LIGHT-POST® devices have been designed to support and secure the coronal restoration, in case of insufficient residual tooth structure (<4 mm).

D.T. LIGHT-POST® drills are intended to prepare the tooth for insertion of a post that will support and secure the coronal restoration, where insufficient residual tooth structure (<4 mm) exists.

These products should only be used in hospitals, clinics or dental office with qualified dental personnel (dentists).

INDICATIONS FOR USE

D.T. LIGHT-POST® devices are indicated in case of insufficient residual tooth structure (<4mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

D.T. LIGHT-POST® devices are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of **D.T. LIGHT-POST®** devices for pregnant women, breastfeeding women and children under 18 years old (for posts only). It is not recommended for such use.

If the device is used in a way different than the one recommended in the DFU, loosening, breakage of the post or extraction of the tooth might occur.

D.T. LIGHT-POST® posts:

- **D.T. LIGHT-POST®** post is single use device. It must not be re-used to avoid the risk of contamination.
- **D.T. LIGHT-POST®** posts can be cleaned with alcohol. They can neither be sterilized nor disinfected in thermodisinfectant.
- The post must be sterilized individually before insertion into the canal the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g EN ISO 11607-1),
 - ✓ Autoclave: type B complying with regulation standard in force (e.g EN 13060),
 - ✓ Sterilization temperature: 134°C – sterilization time: 18 minutes,
 - ✓ 1 cycle only.
- The post must be immersed in isopropanol for 1 min then dried gently with air syringe for 15s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- The use of rubber dam is recommended.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.
- For profoundly flared or ovoid canals, the use of accessory post is indicated for best results.
- A **D.T. LIGHT-POST® posts** re-access kit is available if needed. Contact your RTD distributor.

D.T. LIGHT-POST® drills:

- **D.T. LIGHT-POST®** drills are supplied non-sterile. Disinfect and sterilize **D.T. LIGHT-POST®** drills before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfectant.
- Inspect **D.T. LIGHT-POST®** drills visually before each use: change to new drills when they are damaged, corroded or dull.
- The UNIVERSAL drill may be used only to remove root canal filling partially, and the FINISHING drills may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the **D.T. LIGHT-POST®** drills stated is 15 uses, following this the dental surgeon is informed that the drill should be replaced.
- **D.T. LIGHT-POST®** drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,
 - ✓ Autoclave: type B complying with regulation standard in force,
 - ✓ Sterilization temperature: 134°C – Sterilization: 18 minutes,
 - ✓ Store the sterilized components in a dry, dust-free place,
 - ✓ If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

ADVERSE REACTIONS

- Debonding of the post
- Root fracture
- Fracture of the post

STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
2. Select the correct **D.T. LIGHT-POST®** post size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the **D.T. LIGHT-POST® UNIVERSAL** drill selected (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the **D.T. LIGHT-POST® FINISHING** drill corresponding to the selected **D.T. LIGHT-POST®** post size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal
6. For **D.T. LIGHT-POST®** post, position the colored O-ring as a cutting mark and then remove post. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post. Remove the O-ring.
7. Immersed or soaked the post in isopropanol for 1 minute and then gently air dry for 15 s at least.
8. Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.
9. Apply 2 coats* of the primer with a thin brush (e.g. Composibrush RTD) in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
11. a. If using a dual cure resin cement for both cementation and core build up such as **CORECEM™** (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post

immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic form matrix, or directly on the tooth.

11. b. If using separate cement and a core build up resins: Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
- 11c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
- 11c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the post be covered by core composite.

STORAGE

Store at ambient temperature. Check the packaging of the product before using it. In case of damage, do not use the product.

DISPOSAL

Dispose in a safe manner in accordance with local/national regulations.

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TERMS AND SYMBOLS IDENTIFICATION

	Manufacturer
	Do not re-use
	Batch code
	Catalogue number
	Consult Instructions for use

	Caution
	Number of unit / Content

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INSTRUCTIONS FOR USE

DESCRIPTION

D.T. LIGHT-POST ILLUSION® X-RO® are dental posts used for reconstruction of devitalized teeth. They are made of X-RO yellow translucent fibers bound in a colored resin Epoxy-vinyl ester matrix, which becomes translucent in the mouth. The post color will re-appear on command under cold air, water spray if re-treatment is required.

The **D.T. LIGHT-POST ILLUSION® X-RO®** is a smooth-surfaced post, apically tapered with a length of 20 mm. The optimized cylindrico-conical endo-design limits the removal of dentin, thereby increasing the strength and clinical quality of the restoration.

The radiopacity level of the **D.T. LIGHT-POST ILLUSION® X-RO®** is 340 % Al equivalent (ISO 4049).

The devices are available in 4 sizes and are used in association with **D.T. LIGHT-POST®** drills. A universal pre-shaping drill to perform a first hole and a finishing drill associated to each post size according to an O-ring color part (See Figure 1).

D.T. LIGHT-POST® drills are associated with:

- **D.T. LIGHT-POST** posts
- **D.T. LIGHT-POST ILLUSION X-RO** posts


Figure 1 : **D.T. LIGHT-POST ILLUSION® X-RO®**.

COMPOSITION

D.T. LIGHT-POST ILLUSION® X-RO® posts are fiber-reinforced polymer, made of XRO fiber (80%) and epoxy resin matrix (20%) which contains colored pigment.

D.T. LIGHT-POST® drills are made of stainless steel

INTENDED USE

D.T. LIGHT-POST ILLUSION® X-RO® posts has been designed to support and secure the coronal restoration in case of insufficient residual tooth structure (<4 mm).

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INDICATIONS FOR USE

D.T. LIGHT-POST ILLUSION® X-RO® posts and associated drills are indicated in case of insufficient residual tooth structure (<4mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

D.T. LIGHT-POST ILLUSION® X-RO® posts and associated drills are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of **D.T. LIGHT-POST® ILLUSION® X-RO®** posts and associated drills for pregnant women, breastfeeding women and children under 18 years old (for posts only). It is not recommended for such use.

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D.T. LIGHT-POST® ILLUSION® X-RO® posts:

- The post must be immersed or soaked in isopropanol for 1 minute and dried gently with air syringe for 15 s at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.
The post can neither be sterilized nor disinfected in thermoisinfector and the use of autoclave is not recommended.
- Posts are single use devices. They must not be re-used to avoid the risk of contamination.
- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
- The use of rubber dam is recommended and can momentarily prevent the post to change color.
- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.

D.T. LIGHT-POST® drills:

- **D.T. LIGHT-POST®** drills are supplied non-sterile. Disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfector.
- Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull.
- The universal drill may be used only to remove root canal filling partially, and the finishing drills may be used only to prepare root canal.
- The use of rubber dam is recommended.
- The lifetime of use for the **D.T. LIGHT-POST®** drills stated is **15 uses**, following this the dental surgeon is informed

- that the drill should be replaced.
- Drills must be sterilized individually in an autoclave with the following settings**:
 - ✓ Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use,
 - ✓ Autoclave: type B complying with regulation standard in force,
 - ✓ Sterilization temperature: 134°C - sterilization time: 18 minutes,
 - ✓ Store the sterilized components in a dry, dust-free place,
 - ✓ If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

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STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
2. Select the correct **D.T. LIGHT-POST® ILLUSION® X-RO®** size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the **D.T. Preshaping Drill** (rotation speed 1000-2000 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the **D.T. Finishing Drill** corresponding to the selected **D.T. LIGHT-POST® ILLUSION® X-RO® post** size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal
6. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
7. Clean the post with alcohol.
8. **Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.**
9. Apply 2 coats* of the primer with a thin brush (e.g. Composibrush RTD) in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
11. a. If using a dual cure resin cement for both cementation and core build up such as **Corecem™** (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic form matrix, or directly on the tooth.
11. b. If using separate cement and a core build up resins : Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).
11. c.i. If using a self-etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
11. c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the

post be covered by core composite. In cases where a prosthetic crown is not indicated, a minimum of 1 mm of composite restorative should cover the post and that should be an appropriately-shaded, polishable, aesthetic composite. If the **ILLUSION** post cannot be adequately covered with composite, then a non-**ILLUSION** version can be utilized.

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