Calibra® Ceram

Adhesive Resin Cement

DIRECTIONS FOR USE	ENGLISH	2
Résine composite de collage MODE D'EMPLOI	FRANÇAIS	4
Cemento composito adesivo ISTRUZIONI PER L'USO	ITALIANO	6
Adhäsiver Kompositzement GEBRAUCHSANWEISUNG	DEUTSCH	8
Cemento adhesivo de resina INSTRUCCIONES DE USO	ESPAÑOL	10
Cimento resinoso adesivo INSTRUÇÕES DE UTILIZAÇÃO	PORTUGUÊS	12
Adhesief kunststofcement GEBRUIKSAANWIJZING	NEDERLANDS	14
Adhezinis dervinis cementas NAUDOJIMO INSTRUKCIJOS	LIETUVIŲ K.	17
Adhezīvs sveķu cements LIETOŠANAS INSTRUKCIJA	LATVIEŠU	19
Cimento Resinoso Adesivo INSTRUÇÕES DE USO	PORTUGUÊS (BRASIL)	21

EC REP
Dentsply DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz, Germany
Tel.: 49-7531-583-0
www.dentsplysirona.com

Calibra® Ceram Adhesive Resin Cement

DIRECTIONS FOR USE – ENGLISH

For dental use only. USA: Rx only.

1. PRODUCT DESCRIPTION

Calibra Ceram Adhesive Resin Cement is a visible light cured, dual cured or self cured high strength resin cement which contains fluoride. Calibra Ceram Cement is designed to be used following the application of a compatible dentin/enamel adhesive. Although compatible with numerous indirect restorative materials, there are special instructions for use with light-transmissible restorations.

1.1 Indications

Cementation of indirect restoratives including: Ceramic, composite and metal-based inlays, onlays, crowns, bridges, and posts.

1.2 Contraindications

- Calibra Ceram Cement is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins or any of the components.
- Calibra Ceram Cement is contraindicated for direct application to dental pulp tissue (direct pulp capping).

1.3 Delivery forms (Some delivery forms may not be available in all countries)

Calibra Ceram Cement is available in:

- a convenient 4.5g dual-barrelled syringe
- 5 shades: translucent, light, medium, bleach, and opaque

1.4 Composition

Urethane Dimethacrylate; Di- and Tri-Methacrylate resins; Phosphoric acid modified acrylate resin; Barium Boron FluoroAluminosilicate Glass; Organic Peroxide Initiator; Camphorquinone (CQ) Photoinitiator; Phenylone Oxide Photoinitiator; Accelerators; Butylated Hydroxy Toluene; UV Stabilizer; Titanium Dioxide; Iron Oxide; Hydrophobic Amorphous Silicon Dioxide. Particles of inorganic filler range from 16nm to 7µm, average particle size 3.8µm, total filler 46.3% by volume.

1.5 Compatible adhesives

Calibra Ceram Cement is chemically compatible with conventional methacrylate-based dentin/enamel adhesives including Dentsply Sirona adhesives designed for use with dual-cured resin based materials for all applications. There are special adhesive application instructions for use with light transmissible restorations. See complete directions for use of selected adhesive and step-by-step instructions below. The use of other dentin and enamel adhesive systems is at the discretion and sole responsibility of the dental practitioner.

1.6 Compatible Restoration Substrates

Following application of adhesive as directed (see Directions for Use of selected adhesive product) Calibra Ceram Cement may be used with low and high strength Ceramic materials of construction, e.g., etchable, silica based feldspathic and lithium disilicates, non-etchable, non-silica based alumina and zirconia, with special instructions for light-transmissible restorations ≤ 2.5mm thickness. It is also suitable for use with all resin-based composite, and base, noble and high-noble metallic alloys, full cast or porcelain fused to metal. Calibra Ceram has a chemical affinity for these metals.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other sections of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Calibra Ceram Cement is acidic in nature and contains polymerizable acrylate and methacrylate monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
 - Avoid Eye Contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - Avoid Skin Contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and alcohol and wash well thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - Avoid Contact with Oral Soft Tissues/Mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
 - Avoid Contact with Oral Soft Tissues/Mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
2. Use of the Light-Transmissible curing technique with restorations that are greater than 2.5mm in any thickness, opaque materials or materials that impede light transmission may result in low bond strength, post-operative sensitivity or premature bond/restoration failure. See Step-by-Step instructions below and select compatible adhesive and cementation technique.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Insufficient data exist to support the use of Calibra Ceram Adhesive Resin Cement as a core build up, filling material, base or liner.
3. Insufficient data exist to support the use of Calibra Ceram Adhesive Resin Cement for veneer cementation. Use of Calibra® or Calibra® Veneer Esthetic Resin Cement is recommended.
4. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
5. Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
6. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
7. The syringe cannot be reprocessed. To prevent the syringe and from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringe is handled with clean/disinfected gloves. As additional precautionary measure, syringe may be protected from gross debris but not from all contamination by applying a protective barrier.
8. Calibra Ceram Cement behaves differently intraorally than in ambient operatory conditions. The set of Calibra Ceram Cement is accelerated by the warmth and moisture of the oral environment and/or

ambient or operatory light. After placing Calibra Ceram Cement in contact with tooth structure, e.g., within endodontic post space or in inlay/onlay preparations, immediately seat restoration. Any delay may allow polymerization to begin, which may prevent complete seating of the restoration.

9. Cement at the margins may appear set before cement under the restoration is set. Do not move, torque or disturb restoration until final set of the cement (5 minutes from the start of mix or in the case of light-transmissible restorations completion of light curing).
10. Exercise caution when placing multiple units not to exceed available working time. Load and insure complete seating of each unit prior to loading successive units. The use of a fresh mix of cement is recommended if available working time is approached (see Interactions).
11. Containers should be tightly closed immediately after use. Syringe should be tightly closed by replacing the original cap immediately after use.
12. Calibra Ceram Cement should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause syringe rupture.
13. The Calibra Ceram Cement is a light-curable material. Proceed immediately once materials have been dispensed or protect from ambient light.
14. Use only in well ventilated areas.

Interactions

- Eugenol containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
- Contact with some astringent solutions may interfere with hardening of the polymeric components of the material.
- As with any dual cure resin cement system, the use of an adhesive system can shorten working time. This effect should be investigated in the laboratory prior to clinical use.
- Variable in-vitro data exist regarding use of Calibra Ceram Adhesive Resin Cement in dual cured or self-cured (limited or no light curing) applications in conjunction with some light-cured-only adhesives. Chemical/Product incompatibility may adversely affect product efficacy, leading to premature restoration failure.
- Calibra Ceram contains an optimized photoinitiator combination. Use of high power, broad spectrum output curing lights during cleanup may produce unexpected results, rapid hardening and shortened cleanup times. Monowave output LED lights with a single peak output around 470nm are recommended. Check curing light effect on mixed cement in the laboratory prior to clinical use.

2.3 Adverse reactions

1. Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous membranes: inflammation (see Warnings).
2. Product may cause pulpal effects. (See Contraindications)
3. Residual hardened excess cement may lead to soft tissue injury or irritation (See Interactions and Step-by-Step instructions).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Keep out of direct sunlight and store in a well ventilated place at temperatures between 2°-24°C/35°-75°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS

The following preparatory steps apply to cementation of all types of indirect restorations:

1. Following removal of the temporary restoration and any remaining temporary cement, clean enamel and dentin as directed by the adhesive manufacturer's directions.
2. Check the fit and esthetics of the restoration. **Technique Tip:** Occlusal adjustment of inlays/onlays is best accomplished after final cementation.
3. The final shade of cured material is color stable. For color identification, please refer to the Try-In Paste section below.

TRY-IN PASTE (OPTIONAL)

- Dispense the appropriate shade of Calibra® Try-In Paste (available separately) from the syringe onto a clean mixing pad or glass slab. Load paste onto internal surfaces of the restoration and gently seat onto preparation. Clean excess with a cotton pellet and/or blunt explorer. Shades may be blended to achieve optimum esthetics. **Important Technique Tip:** Factors such as restoration shape, texture, thickness and shade(s) of the underlying and surrounding tooth structure contribute to final perceived shade. Try-In Paste is a guide for cement shade range selection only. **NOTE:** The try-in paste will not polymerize, thus offers unlimited work time.
- Once fit and esthetics are verified, thoroughly rinse Try-In Paste from restoration and preparation surfaces using water.

TREATMENT OF THE RESTORATION

Metal Restorations

- Internal surfaces of restorations should be clean and dry prior to cementation. Internal surface microetching (sandblasting with 50µ alumina) of metal surfaces of the restoration is recommended.
- **Technique Tips For Maryland Bridges:** The use of 180° wrap of wings, rest seats, parallelism and slots or grooves is necessary for appropriate retention. The metal wings of the Maryland Bridge should be perforated, electrolytically etched, laboratory chemically etched, or mechanically sandblasted with 50µ alumina.

Ceramic/Composite Restorations

- Follow the dental laboratory or restoration manufacturer's instructions for pre-treatment, if required. Restorations designed to be silanated or if the internal silanated surface has been disturbed during try-in, apply Calibra® Silane Coupling Agent (available separately) as outlined in the manufacturer's Directions for Use.

TOOTH CONDITIONING/DENTIN PRETREATMENT/ADHESIVE APPLICATION

- Following restoration try-in, rinse preparation thoroughly with water spray and air dry.
- Proceed with adhesive application as directed in the complete directions for use.

Delayed or Reduced Adhesion due to adhesive incompatibility

- Select compatible light cure or dual cure adhesive

- Always refer to and follow adhesive manufacturer's complete directions for use for selection and application of the appropriate light cure or dual cure adhesive. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately following adhesive application to cementation with Calibra Ceram Cement, dual or self-cured dual barrel mix as indicated below.
- Calibra Ceram Cement has been formulated to be compatible with specific* light cured Dentsply Sirona Adhesives without the need for addition of Self Cure Activator (see Complete Directions For Use). Some products may not be available in all countries.

*Prime&Bond elect® Universal Dental Adhesive; Prime&Bond active® Universal Adhesive; Prime&Bond universal® Adhesive

Delayed or Reduced Adhesion due to adhesive incompatibility

- Only Prime&Bond elect Adhesive, Prime&Bond active Adhesive or Prime&Bond universal Adhesive are indicated for use without addition of Self Cure Activator.
- Insufficient data exist to support use of other Dentsply Sirona adhesives without Self Cure Activator. Follow selected adhesive's Dual Cured Cementation procedures.

• Follow adhesive application, air drying (solvent evaporation) and light curing instructions of selected adhesive prior to restoration cementation with Calibra Ceram Cement for all types of restorations.

Special Instructions for Use with Light Transmissible restorations, including most* in-office CAD-CAM milled restorations:

- Final restoration must be light transmissible, ≤ 2.5mm thickness, with light curing of Calibra Ceram Cement through the restoration after seating and excess cement cleanup.
- Select and apply Prime&Bond elect* Adhesive, Prime&Bond active* Adhesive or Prime&Bond universal* Adhesive. Light curing of applied adhesive may be accomplished before or after seating restoration with Calibra Ceram Cement.
- Following application and air drying (solvent evaporation), proceed immediately to cementation, dual cure cleanup and Light Transmissible Visible light curing through the restoration.

*Zirconium-based restorations are generally not light-transmissible. Cure adhesive on the preparation prior to cementation.

Technique Tip: The adjacent teeth and/or the external surfaces of the restoration may be lubricated with a water soluble medium to ease clean up of excess cement.

3.1 Crown/bridge/inlay/onlay/post cementation

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/ Dentin Pretreatment/Adhesive Application. Refer to manufacturer's Directions for Use for selected adhesive. Adhesives (plus activators, if applicable) compatible with both light-cured, CQ initiated methacrylates and self-cured methacrylates are recommended. Apply adhesive (plus activator if applicable) to the internal bonding surface of the restoration if directed by adhesive manufacturer's directions for use.

3.1.1 Cementation technique

Danger of injury due to excessive force

- Apply slow and steady pressure on the syringe
- Do not use excessive force – dual barrel syringe rupture may result

Dual Barreled syringe dispensed

- Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
- Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip with the v-shape notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.
- Gently depress syringe plungers to begin the flow of material. DO NOT USE EXCESSIVE FORCE. If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.

3.1.2 Apply a uniform layer of cement on the entire internal surface of the restoration directly from the mixing tip. For inlays/onlays, it may be helpful to apply a thin layer of cement to the internal portions of the tooth preparation to avoid any porosity or voids. **Technique Tip:** The mixing tip may be bent slightly to allow direct intraoral access for placement of cement into preparations with internal anatomy. For endodontic post spaces, use of a Lentulo Spiral or metal file to aid placement in the post space is recommended. At room temperature, Calibra Ceram Cement offers a minimum work time of 2 minutes when protected from ambient light. If loaded into restoration lined with adhesive or adhesive/activator mixture, room temperature work time (restoration seating time) will be less (see Precautions).

3.1.3 Immediately seat the restoration in the mouth with gradual pressure. Verify complete seating. A gentle rocking or vibratory motion may be helpful to insure optimal seating.

Incomplete cementation

- Intraoral work time is reduced due to warmth and contact with adhesive
- Assure restoration is fully seated within 1 minute
- Placement of multiple units may require multiple mixes of cement

3.1.4 Protect restoration from contamination and movement until the final set of the cement (5 minutes from start of mix or in the case of light-transmissible restorations, completion of light curing).

Inadequate retention due to insufficient stabilization

- Stabilize restoration during cleanup
- Protect restoration from movement during setting
- Allow cement to self-cure without disturbing for 5 minutes

3.2 Cleaning marginal excess Calibra Ceram Cement

All shades of Calibra Ceram Cement are radio-opaque, with a 1mm radio-opacity equivalent to 2mm radio-opacity of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. Thus 1mm of material having a radio-opacity equivalent to 1mm of aluminum has a radio-opacity equivalent to that of dentin.

3.2.1 Self-cure cleanup

The excess cement will reach the "gelled" state after approximately 1-2 minutes in the mouth, allowing easy removal. Excess cement will remain in the "gelled" state for approximately 1 minute. If exposed to directed operatory light, "gel" state may be reached sooner and remain "gelled" for a shorter period. Immediately after reaching the "gelled" state, floss interproximally, only in the direction of restoration seating, to remove excess cement. Complete excess cement removal using an instrument such as a rubber tip, a scaler or an explorer. **NOTE:** Cement within the crown has not yet set. Do not move, torque, or disturb the crown during cleanup. **Technique Tip:** Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.2.2 Optional dual cure cleanup

Due to the dual cure property of Calibra Ceram dual barrel dispensed cement, the operator has the option of utilizing a curing light to facilitate cleanup. Light curing to facilitate cleanup must be accomplished within the first minute following intraoral insertion. Excess cement cleanup may begin immediately following a brief exposure with the curing light. Conventional powered quartz tungsten halogen or LED lights producing only one peak wavelength around 470nm are recommended. Use of narrow spectrum output lights or high power, broad spectrum output halogen or LED lights may produce unexpected results. Light curing mixed cement after one minute, or continuously at one spot for more than 5 seconds, at any time, will cause cement to adhesively set, making cleanup difficult.

Immediately after seating is verified, briefly light-cure excess cement at the margins by constantly moving the curing light tip around the margins for no more than 5 seconds per surface. Excess cement will reach a "gelled" state after this brief (< 5 second) cure.

Residual excess cement - overcure during cleanup

- Monowave (470nm) LED or conventional halogen lights are recommended
- High power, dual or broad spectrum lights may cause premature hardening of excess cement
- Evaluate curing light/cement interaction in laboratory prior to clinical use

The excess cement will remain in the "gelled" state for approximately 45 seconds following light exposure. All excess cement must be removed before final self-cure set is achieved, as outlined above. **NOTE:** Cement within the crown has not yet gelled or set. Do not move, torque, or disturb the crown during cleanup. Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.3 Curing, finishing and dismissal

Non-Light-Transmissible restorations (Self-Cure)

For metallic, thick or heavily opaqued ceramic or composite, or restorations that otherwise impede the transmission of light, once cleanup is completed and restoration is stabilized, allow Calibra Ceram Cement to self-cure without disturbing for 5 minutes from start of mix. Protect restoration from contamination and movement during the setting time. Following the self-cure set, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

Inadequate retention due to insufficient stabilization

- Stabilize restoration during cleanup
- Protect restoration from movement during setting
- Allow cement to self-cure without disturbing for 5 minutes

Light-Transmissible restorations (Light-Cure)

For most non-metallic, light-transmissible ceramic or composite restorations, Calibra Ceram Cement may be visible light cured. **Important Technique Note:** If following the Special Instructions for In-office Milled Ceramic Restorations, utilizing uncured adhesive technique, Calibra Ceram Cement must be visible light cured through the restoration after seating. Once cleanup is completed and restoration is stabilized, light cure all areas of the restoration using a visible light, curing unit (output spectral peak including 470nm, broad spectrum output lights recommended), with a minimum output of 550mW/cm² for 20 seconds from each direction – buccal, lingual and occlusal. Following the light-curing, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

Inadequate polymerization due to insufficient curing

- Check compatibility of curing light
- Check curing cycle
- Check curing output before each procedure

Important Technique Tips:

- When simultaneously cementing multiple single units or bridgework, it is recommended to employ the light-cure cleanup on one or two adjacent units only, allowing other units' excess to self-cure, providing ample cleanup time.
- Clean excess cement from metal instruments immediately as set cement will adhere to the instrument.

3.4 Finishing and polishing

- Removal of resin flash is best accomplished with the Enhance® Finishing System of points, cups and discs. The Enhance® System will remove flash and finish restoration margins without removal of or trauma to the enamel.
- Polish final restoration using Enhance® PoGo® Polishing system or Prisma® Gloss™ Polishing Paste and Prisma® Gloss® Extra Fine Polishing Pastes. (See complete Directions for Use supplied with polishing product chosen)

4. HYGIENE

Cross-contamination – Infection

- Do not reuse single use products. Dispose of in accordance with local regulations.
- To prevent the syringes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringes are handled with clean/disinfected gloves. Do not reuse syringes if contaminated.
- Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.

For dual-barreled syringe, remove used mixing tip and discard appropriately. Replace original syringe cap prior to storing.

- To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination.

Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard composite material that has been in contact with any fluid or non-sterile instrument.

- Repeated disinfection may damage label.

NOTE: Destruction of the label by vigorous wiping. Wipe syringe gently.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard is used: "YYYY-MM-DD"
2. The following numbers should be quoted in all correspondence:

- Reorder Number
- Lot number
- Expiration date

© 2019 Dentsply Sirona Inc. All Rights Reserved. 507104 (R 9-27-19)

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