Gutta-Percha Points

EN	Instructions For Use	2
FR	Instructions d'utilisation	9
IT	Istruzioni per l'uso	16
DE	Gebrauchsanweisung	23
ES	Instrucciones de uso	30
PT	Instruções de utilização	37
NL	Gebruiksaanwijzing	44
SV	Användarinstruktioner	51
DA	Brugsanvisning	58
NO	Brukerinstruksjoner	65
LT	Naudojimo instrukcija	72
LV	Lietošanas norādījumi	79
ET	Kasutusjuhised	86
CS	Návod k použití	93
SK	Návod na použitie	100
HU	Használati utasítás	107
RO	Instrucțiuni de utilizare	114
SR	Uputstva za upotrebu	121
BS	Uputstva za upotrebu	128
CNR	Uputstva za upotrebu	135
SQ	Udhëzimet për përdorim	142
EL	Οδηγίες χρήσης	149
TR	Kullanım Kılavuzu	156
PT-BR	Instruções de uso	163
KO	사용 지침	170

Gutta-Percha Points

Instructions
For Use

FOR DENTAL USE ONLY
 FOR PROFESSIONAL USE ONLY
 SINGLE USE
 NON-STERILE – TO BE DISINFECTED BEFORE USE

INSTRUCTIONS FOR USE

RECIPROC® MINIMA CONFORM FIT® GUTTA-PERCHA POINTS (REF. MNSRMGPF)

1 PRODUCT DESCRIPTION

- Product description:
Gutta-Percha Points are endodontic obturation devices made of radio-opaque synthetic gutta-percha rubber used to fill the root canal after cleaning and shaping. These devices are single-use and available only in non-sterile version.
- Clinical benefits:
The clinical benefit to the patient, when performing an obturation procedure with Gutta-Percha Points is to restore tooth integrity and masticatory function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
The products are indicated for treatment of endodontic disease.
- Intended purpose:
The products are intended to fill the shaped, cleaned, and irrigated root canal space.
- Intended patient population:
Patients who need root canal treatment.
- Intended user and environment of use:
Endodontic devices are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

Code	Product	Tip size Taper	Packaging and content
MNSRMGPF	RECIPROC® Minima Conform Fit® Gutta-Percha Points	R20 R25 R35 R45 Taper variable	Box of 60 points

Code	Product	Tip size Taper	Packaging and content
MNSRMGPF00AST	RECIPROC® Minima Conform Fit® Gutta-Percha Points	R20 R25 R35 R45 Taper variable	Box of 60 points including: 12 points of size R20 24 points of size R25 12 points of size R35 12 points of size R45

5 COMPOSITION

- Components:
 - Zinc Oxide (45 - 55%)
 - Synthetic Gutta-Percha (20 - 30%)
 - Barium Sulfate (15 - 20%)
 - Titanium Dioxide (1 - 5%)
 - Yellow Colorant (< 2%)
 - Red Colorant (< 2%)
 - Remaining components (< 10%)

Physical and chemical properties

- Radiopacity: > 3 mm of aluminum equivalent in thickness
- Melting Temperature: < 100°C (212°F)
- Nominal length: ≥ 29 mm (1.1 in)

6 COMPATIBLE DEVICES

Gutta-Percha Points are used in combination with:

- Spreaders and Pluggers. These instruments are used to assure the gutta-percha point is distributed throughout the entire root canal system.
- Calamus® and Gutta-Smart™. They are heating devices that are used to extrude gutta-percha material into the root canal system.
- Gutta-Condensor. Gutta-Condensor is a thermal condenser used to soften and fuse gutta-percha points together.
- Endodontic sealers, such as AH Plus® Sealer. AH Plus® Sealer and ProRoot® Bio Sealer are used to coat root canal walls and remains at the interface between the root canal and the gutta-percha material.
- Shaping files and Paper Points. They are devices used to shape and dry the root canal system.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions for Endodontic Products (see section [9]) to minimize the following risks to the device, the patient and/or the user:
 - Cross-contamination. Gutta-Percha Points obturators are intended for single use only. Single-use products shall not be processed and reused. Reusing these products increases the risk of cross-contamination and/or breakage.
 - Use of the devices might lead to the following residual risks:
 - Dissatisfaction of the user and delay in treatment / patient discomfort.
 - Retreatment (initial infection never healed).
 - Serious infection / serious health consequences for the patient (e.g. biological hazards such as infection, toxic or allergic reaction for the patient) have never been reported and are very unlikely but cannot be totally excluded.
 - Fracture of the tooth due to obturation instruments.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Inspect the packaging and product before use: if the packaging or the product is damaged, do not use the product and discard it.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Do not sterilize with heat.
- Do not use disinfecting solutions containing Phenol or any products which are not compatible with gutta-percha.

7.3 Adverse reactions / Undesirable side effects

- Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

- Keep away from sunlight and heat.
- Maintain under room temperature ($\leq 30^{\circ}\text{C}$ (86°F)).
- Store in a dry area.

8 STEP BY STEP INSTRUCTIONS

The described Gutta-Percha Points are part of a complete system and should only be used in combination with the corresponding shaping files and paper points. They are identifiable with their system brand, size name and color code.

The Gutta-Percha Points shall only be used in a clinical or hospital environment, by qualified dental professionals. After shaping and cleaning the root canal, the obturation technique must be set by the professional.

This procedure includes the use of Gutta-Percha Points and sealer.

- 1) Gutta-Percha Points selected shall come from the same shaping system and correspond to color code and shaping file size of the last instrument used during canal preparation. The point shall fit the apical size of the last instrument used and a small friction should be felt when withdraw (tug-back).
- 2) The points are to be disinfected prior to use (see section [9]).
- 3) A heated instrument can be used to soften the material and ensure perfect sealing and dimensional adaptation of the root canal.
- 4) In the use of a heated instrument, follow the IFU of said instrument.
- 5) Use an appropriate tool to cut the Gutta-Percha Points, including color coding ring, if present, after the insertion of the Gutta-Percha Points in the root canal. The colored ring shall not remain in the root canal after insertion.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- Gutta-Percha Points are intended for Single Use only.
- Gutta-Percha Points are non-sterile devices. Please, follow the below instructions for disinfection:
 - Immerse the Gutta-Percha Points in 5.25% NaOCl for 1 to 5 minutes at ambient temperature.
 - Remove all the bubbles from the surface of the Gutta-Percha Points.
 - Gently wipe the Gutta-Percha Points with sterile gauze moistened with 70% isopropyl alcohol.
 - Let the Gutta-Percha Points air dry.
- For other products listed in this IFU, please follow processing instructions from the manufacturer.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Information will be available at both the European database on medical devices (EUDAMED) and the IFU website, where it will be linked to the Basic UDI-DI.

EUDAMED

<https://ec.europa.eu/tools/eudamed>

IFU

dentsplysirona.com/ifu

The Basic UDI-DI number is: ++J00310043DX

15 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB
United Kingdom

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Dentsply Sirona Europe GmbH