Anti-bacterial Retainer Brite® Dental Appliance Cleaner

EN

INSTRUCTIONS FOR USE

INTENDED USE:

Product is an anti-bacterial dental appliance cleaner intended to remove bacteria and clean a variety of removable dental appliances.

INDICATIONS FOR USE:

Product is indicated for patients who are using removable dental appliances such as retainers, aligners, dentures, sports mouthguards, nightguards, bruxing and TMJ appliances, which need to be cleaned at home.

PERFORMANCE CHARACTERISTICS AND CLINICAL BENEFITS:

Retainer Brite® Dental Appliance Cleaner is specially formulated to help protect removable appliances by helping remove plaque and tartar, killing common odor-causing bacteria* and preserving clarity and cleanliness of the appliance.

Clinical benefits of orthodontic treatment may include overall improvement of tooth health and correction of malocclusion.

CONTRAINDICATIONS:

This product is contraindicated for patients allergic to Persulfates.

WARNINGS:

Contains persulfates which are a known allergen.

Keep out of reach of children.

Use only as directed. If under the care of a dental professional, follow the instructions provided by your dentist.

If product comes into contact with eyes, wash eyes thoroughly with water for several minutes.

Do not put tablets or cleaning solution in mouth. If ingested, seek medical assistance immediately. In the USA, please call Poison Control Center at 1-800-222-1222. In the UK, please call NHS 111.

Report serious incidents that have occurred in relation to the device to device manufacturer and the Competent Authority in Member State of residence, or other applicable regulatory agency.

PRECAUTIONS:

**Do not allow wire appliances to soak for more than 15 minutes. Remove promptly and rinse thoroughly with water.

Wash hands thoroughly before and after handling tablets.

Do not mix with other cleaners.

DO NOT REUSE, once dissolved tablet cannot be reused. Solution will not clean as intended if reused.

ADVERSE REACTIONS:

Allergic reactions.

STEP BY STEP INSTRUCTIONS:

Use daily for best results.

STORAGE:

Keep Dry

SAFE DISPOSAL OF THE DEVICE:

Dispose of down the drain, as you would your mouthwash or toothpaste. Flush drain with plenty of running water.

*Kills 99% of common odor-causing bacteria including: *P. aeruginosa*, *S. aureus*, *E.coli*, *F. nucleatum* (ASTM E2315-03), *P. aeruginosa*, *S. aureus*, *E.hirae*, *E. facium* (EN 13727 and EN 14561). Independent test lab data on file.

INGREDIENTS: Sodium Bicarbonate, Citric Acid, Sodium Carbonate, Potassium Persulfate Compound, Corn Syrup Solids, Sodium Percarbonate, Sodium Sulfate (may contain silica), Sorbitol, Tetraacetylenediamine (TAED), PEG-180, Sodium Lauryl Sulfoacetate, Flavor, PEG-8, Magnesium Stearate, FD&C Blue #1 Aluminum Lake (CI 42090), FD&C Blue #2 (CI 73015)

Manufacturer:
Dentsply Sirona Orthodontics Inc.
7290 26th Court East
Sarasota, FL 34243 USA

Authorized EU Representative
Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Authorized Switzerland Representative
Maillefer Instruments
Holding Sàrl
Chemin du verger 3
CH-1338 Ballaigues

www.dentsplysirona.com/ifu

MKT-DFU0048-ML Rev. 7 Issued 01-2024
© 2024 Dentsply Sirona

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